EN

DA

DE

EL

ES

FI

FR

IT

JA

NL

NO

PL

PT

SV

TR

INSTRUCTIONS FOR USE

The logo for 'isodry' features a stylized grey 'i' with a curved line above it, followed by the word 'isodry' in a bold, lowercase sans-serif font. A registered trademark symbol (®) is located to the upper right of the 'y'.

Dental Isolation System

Better Isolation = Better Dentistry®

TABLE OF CONTENTS
INDHOLDSFORTEGNELSE
INHALTSVERZEICHNIS
ΠΙΝΑΚΑΣ ΠΕΡΙΕΧΟΜΕΝΩΝ
ÍNDICE
SISÄLLYSLUETTELO
TABLE DES MATIÈRES
INDICE
目次
INHOUDSOPGAVE
INNHOLDSFORTEGNELSE
SPIS TREŚCI
ÍNDICE
INNEHÅLLSFÖRTECKNING
İÇİNDEKİLER

The logo for Isodry, featuring a stylized 'i' with a grey arc above it, followed by the word 'isodry' in a bold, lowercase sans-serif font, and a registered trademark symbol (®) to the right.

Dental Isolation System

Better Isolation = Better Dentistry®

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Isodry® Dental Isolation System

[REF] EID0101 (USA) [REF] EID0103 (EU) [REF] EID0104 (AUS)

This system is intended for use with the single-use, sterile or non-sterile Mouthpieces developed by Isolite® Systems. These Mouthpieces are available in a wide-range of sizes for effective isolation.

Each Isodry Dental Isolation System is comprised of three main components: 1) Titanium Control Base, 2) Vacuum Pipe, 3) Vacuum Hose.

Titanium Control Base

[REF] PID0101

The Titanium Control Base includes two Vacuum Levers (P/N: PIL0205) to adjust suction in the upper and lower areas of the Mouthpiece.

Vacuum Pipe

[REF] PID0103

The Vacuum Pipe channels suction to the Mouthpiece. The Mouthpiece attaches onto the Vacuum Pipe.

Vacuum Hose

[REF] PID0102 (USA)
(1.82 meters / 72 inches)
[REF] PID0104 (EU)
(1.49 meters / 59 inches)
[REF] PID0102 (NZ/AUS)
(1.82 meters / 72 inches)

The Vacuum Hose delivers vacuum suction from the vacuum canister to the Control

Head. The Vacuum Hose connects to the vacuum canister.

Suction Line Plug

[REF] PIL0206

O-ring Maintenance Kit

[REF] AIL0206 (USA/NZ/AUS)
[REF] AIL0223 (EU)

3/8" Y-adapter Kit

[REF] AIL0207 (USA/AUS)

9.5 mm Y-adapter Kit

[REF] AIL0222 (EU/AUS)

Double End Tube Brush

[REF] AIL0204

Mouthpiece Lubricant

[REF] AIL0202

Barrier Sleeves

[REF] AIL0203 (non-sterile)
Sold Separately

Non-sterile Control Head Barrier Sleeves are available for additional protection

from cross-contamination. Use of Barrier Sleeves does not eliminate the need to follow disinfecting and sterilization recommendations for the Control Head and Vacuum Pipe.

Mouthpieces

[REF] CILXX0X (non-sterile)
[REF] CISXX0X (sterile)
Sold Separately

Sterile and non-sterile Mouthpieces are single use only and available in various sizes. The Mouthpieces are constructed of a soft, flexible polymer and are not made with natural rubber latex. They do not contain phthalates or BPA's. For more information on sizes, go to isolitesystems.com.

Intended Use

Isodry is a dental isolation system for clinical use by dental professionals. The system is designed to isolate two quadrants at once, simultaneously providing tongue and cheek retraction and continuous suction. The retraction and continuous suction are the important functions that control the oral environment, reducing contamination and improving bond strength while protecting the airway and soft tissue.

Contraindications

No known contraindications. Do not use without reviewing how to select the correct Mouthpiece size for each patient.

Do not use if unable to position the Mouthpiece correctly.

Side Effects

The use of the Mouthpiece may trigger a pharyngeal reflex (laryngeal spasm) in hypersensitive patients, or if it is not placed correctly within the oral cavity.

How Supplied

Isodry

The Isodry Dental Isolation System contains: Titanium Control Base, 6 Vacuum Pipes, Vacuum Hose, and Instructions for Use.

Accessories: Suction Line Plug, Mouthpiece Lubricant, Y-adaptor Kit, Double End Tube Brush, and O-ring Maintenance Kit.

Mouthpieces

Mouthpieces (sold separately) are supplied in convenient operatory boxes. Each Mouthpiece is individually packaged. Non-sterile Mouthpieces are packaged inside a polyethylene bag while Sterile Mouthpiece are packaged inside a thermoform tray with a Tyvek lid. Sterile Mouthpieces are sterilized using e-beam radiation. When sterilized, yellow indicator dots on both the Tyvek lid and the box turn red. Do not use the Sterile Mouthpiece if the indicator dots are not red.

Installation

Isodry is designed to install easily in most dental operatories (Figure 1). Installation and usage instructions, including a video demonstration, are available at isolitesystems.com.

Step 1 Connect the Vacuum Hose to the vacuum canister by fitting the end of the Vacuum Hose onto an auxiliary HVE port on the vacuum canister. If an auxiliary HVE port is not available, use the provided Y-adaptor, so that the standard HVE hose may remain connected.


Figure 1. Installation

NOTE: For countries outside the USA, installation may need to be completed by a qualified service technician. Please consult your in-country representative for additional information.

Step 2 Connect the Control Head to the head of the Vacuum Hose. Once connected, test the dual Vacuum Levers to confirm suction functionality.

Step 3 Place the Isodry in an open slot on the toolbar when not in use.

How to Operate

Using the Dual Vacuum Levers

Dual Vacuum Levers adjust suction intensity in both the upper and lower quadrants. To adjust, rotate the Left and Right Vacuum Levers to the desired suction intensity for each quadrant. To shut off suction, rotate the Levers to the most downward position.


Figure 2. The Finger Method

Selecting the Mouthpiece

To select the correct Mouthpiece size, measure the patient's interincisal opening. If between sizes, you may select the smaller size.

| Mouthpiece Size | Approx. Interincisal Opening |
|-----------------------------|------------------------------|
| Pediatric (PED) | < 30 mm |
| Extra Small (XSM) | < 30 mm |
| Small (SML) | 30 mm |
| Medium (MED) | 30-45 mm |
| Medium Deep Vestibule (MDV) | 30-45 mm |
| Large (LRG) | > 45 mm |

Table 1. Mouthpiece Sizes

The Finger Method

The finger method may also be used to select the appropriate Mouthpiece size (Figure 2). For additional assistance, watch a video demonstration online at isolitesystems.com.

Attaching the Sterile Mouthpiece

Follow the established protocol in your practice for handling sterile product. Open the thermoform tray by removing the Tyvek lid and place the Mouthpiece in a sterile field. Slide Mouthpiece onto the Vacuum Pipe until the two locking bars engage. Sterile water may be used to lubricate the inside surface of the Mouthpiece sleeve to ease attachment.


Figure 3. Mouthpiece

Attaching the Non-sterile Mouthpiece

Slide the Mouthpiece onto the Vacuum Pipe until the two locking bars engage. Water or mouthpiece lubricant may be used to lubricate the inside surface of the Mouthpiece sleeve to ease attachment.

Placing the Mouthpiece

Prior to placing the Mouthpiece within the patient's oral cavity, moisten the Mouthpiece and the patient's lips.

Step 1

Hold the sides of the Control Head with your index finger and thumb. With your other hand, fold the Cheek Shield onto the Tongue Retractor. Gently slide the folded mouthpiece into the buccal vestibule on the side to receive treatment. Angle the mouthpiece to allow the lower edge of the Tongue Retractor

to move along the buccal edge of the teeth.

Step 2

Move the bite block onto the occlusal surface of the teeth, just distal to the mandibular cuspid. Instruct the patient to “rest gently” on the bite block to secure.

Step 3

Place the Cheek Shield into the buccal vestibule. Move the Isthmus onto the retromolar pad behind the maxillary tuberosity. Adjust the Tongue Retractor in the lingual vestibule as


Figure 4, 5, & 6. Mouthpiece Placement

needed. Move the bite block distally to provide more vertical working room.

Disassembly of the Control Head

Step 1

Press the release button (Figure 7) on the underside of the Titanium Control Base to release and disconnect the Vacuum Pipe.

Step 2

Holding the head of the Vacuum Hose, pull to remove the Titanium Control Base. Do not twist.


Figure 7. Control Head

Disinfection, Cleaning, and Sterilization

Surface Disinfection of the Titanium Control Base

Isolite Systems recommends performing a surface disinfection of the Titanium Control Base between each patient. To disinfect the Titanium Control Base, first remove the Vacuum Pipe.

If the Titanium Control Base is visibly soiled with biomatter or blood, use an intermediate level disinfectant wipe with a claim to inactivate *Mycobacterium tuberculosis*. Otherwise, use a hospital grade disinfectant and follow the manufacturer's instructions for contact time.

Surface Disinfection of the Vacuum Hose

Isolite Systems recommends performing a surface disinfection of the Vacuum Hose between each patient.

Use a hospital grade disinfectant per manufacturer's instructions. If the components are visibly soiled with biomatter or blood, use an intermediate level disinfectant wipe with a claim to inactivate *Mycobacterium tuberculosis*. Barrier Sleeves may also be used and changed between patients.

Cleaning and Sterilization of the Vacuum Pipe

Isolite Systems recommends sterilization of the Vacuum Light Pipe between each patient.

Remove biomatter and other debris from the suction channels using the provided Double End Tube Brush with a neutral pH (6.5-7.5) enzymatic cleaning solution. Do not use solution to hold parts. After cleaning, rinse with tap water, dry, and place components in an autoclave pouch to be sterilized in a steam autoclave. Use validated sterilization cycles (Table 4).

Cleaning and Sterilization of the Titanium Control Base

Isolite Systems recommends following normal infection control protocol for dental devices. Users may sterilize the Titanium Control Base at their discretion. The Titanium Control Base should be sterilized when used with the Sterile Mouthpiece.

Remove Vacuum Levers from the Titanium Control Base. Remove biomatter and other debris from the suction channels using the provided Double End Tube Brush with a neutral pH (6.5-7.5) enzymatic cleaning solution. Do not hold parts in cleaning solution or use ultrasonic cleaners. After cleaning, rinse with tap water, dry, and place components in an autoclave pouch to be sterilized in a steam autoclave. Use validated sterilization cycles (Table 4).

After sterilization, lubricate the Vacuum Lever O-rings with the Dow 111 O-ring lubricant provided in the O-ring Maintenance kit.

Validated Sterilization Cycles

| Cycle | Temperature | Min. Exposure Time | Min. Drying Time |
|-----------|-----------------|--------------------|------------------|
| Standard | 121 °C / 250 °F | 30 min | 30 min |
| Standard | 132 °C / 270 °F | 15 min | 15 min |
| Prevacuum | 132 °C / 270 °F | 4 min | 20 min |
| Prevacuum | 134 °C / 273 °F | 3 min | 20 min |

Table 2. Validated Sterilization Cycles

Maintenance

Titanium Control Base

The Vacuum Lever O-rings require weekly lubrication to operate at peak efficiency. Remove the Titanium Control Base from the Vacuum Hose. To remove the Vacuum Levers, push the levers to the “off” position and keep pushing gently until each lever pops out, or gently pry them out. Coat the O-ring by applying a small amount of the Dow 111 O-ring Lubricant provided in the O-ring Maintenance Kit. Apply an additional amount to the Titanium Control Base retainer rings as well. Re-attach the Vacuum Levers to the Titanium Control Base (Figure 8 and 9).

Vacuum Hose

Clean the inside of the Vacuum Hose regularly according to frequency of use to prevent the accumulation of biofilm. Start by vacuuming water through the Vacuum Hose. Then vacuum a dental hose cleaner through the hose per manufacturer’s instructions. The vacuum hose cleaner should be non-foaming, deodorizing, anti-microbial, and able to dissolve organic debris. Most enzymatic cleaners with pH between 6 and 8 are acceptable. Do not use disinfectants as a vacuum hose

Figure 8. Control Base

Figure 9. Vacuum Levers

Figure 10. Vacuum Hose

cleaner. If further disinfection is necessary, consult your vacuum manufacturer for recommendations.

To lubricate the Vacuum Hose O-ring, remove the Titanium Control Base from the Vacuum Hose. The black O-ring is located approximately 13 mm (½ inch) down the inside of the blue fitting at the head of the Power Hose. (See Figure 10) Swab the O-ring with the Dow 111 O-ring Lubricant provided in the O-ring Maintenance Kit.

Warning

⊗ **WARNING:** Mouthpieces are single use only. Reuse can result in the transmission of infectious pathogens. Mouthpieces must be discarded after use. Mouthpieces cannot be safely resterilized after use.

⊗ **WARNING:** Do not use if the individual Mouthpiece package has been previously damaged or opened. Non-sterile Mouthpieces are individually packaged in polyethylene bags. Sterile Mouthpieces are individually packaged in a thermoform tray with a sealed Tyvek lid.

⚠ **WARNING:** Be sure to use established protocols and patient monitoring when using the Mouthpiece in conjunction with sedation dentistry. Patient management by a certified anesthesiologist is strongly recommended.

Precautions

⚠ **CAUTION:** Do not use the Sterile Mouthpiece if the sterilization indicator dot on the Tyvek lid is not red.

⚠ **CAUTION:** When performing dental procedures requiring high levels of blood evacuation, it is recommended that the Mouthpiece be intermittently flushed with water to prevent coagulation in the Mouthpiece to avoid loss of suction and light transmission.

⚠ **CAUTION:** Do not autoclave Titanium Control Base at temperatures above 135 °C; higher temperatures could damage components and void your product warranty.

⚠ **CAUTION:** Do not hold Titanium Control Base in cleaning solution or use ultrasonic cleaners.

⚠ **CAUTION:** Failure to maintain properly lubricated O-rings can result in damage to the Isodry that is not covered by product warranty.

Guide to Symbols

Limited 1-Year Warranty

Isolite® Systems warrants this product to be free from defects in material and workmanship for a period of one year from the date of original purchase. If the product exhibits such a defect, Isolite Systems will, at its option, repair or replace it without costs for parts and labor. The product must be returned by the customer to Isolite Systems in accordance with current Return Material Authorization procedures. Proof of date of original purchase may be required. This warranty does not cover finishes or normal wear, nor does it cover damage resulting from accident, misuse, dirt, tampering, unreasonable use, failure to provide reasonable and necessary maintenance, service performed or attempted by unauthorized individuals or service agencies, or devices that have been otherwise modified. All implied warranties, including any implied warranty or merchantability or fitness for any particular purpose, are limited in duration to one year from the date of receipt of the product. In no event will Isolite Systems be responsible for consequential damages resulting from the use of this product.

Instructions for Use may also be found at: isolitesystems.com.

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All rights reserved.
Innerlite, Inc.