


Attention!!!

Solutions for ADEC Cascade Owners


Isolite works best when you have two available HVE ports. However, ADEC Cascade units are pre-installed with the following component:



NOTE: It is NOT necessary to replace this component before using Isolite. Call Isolite, 1-800-560-6066 for information.

Your supplier has a replacement option that will make Isolite easier to use.

A second high-speed vacuum will allow you to suction larger pieces of debris not aspirated by the Isolite Mouthpiece.



Contact your local technician to discuss replacement options for your Solids Collector. See graphic for ADEC part numbers.